PlumePen® Pro

Surgical Smoke Evacuation Pencil

Intended Use: PlumePen® Pro is designed for general electrosurgical applications, including cutting and coagulation, and for removing smoke generated by electrosurgery when used in conjunction with an effective smoke evacuation system. The pencil enables the operator to remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the operative site for the desired surgical effect.

Indications for Use:

- To remove smoke plume from the surgical site.
- To remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the target tissue for the desired surgical effect.

Figure 1. Open product packaging and remove PlumePen® Pro assembly. Connect the tubing provided with PlumePen® Pro to a surgical smoke evacuator or other appropriate vacuum source; example shown, (Figure 1). For connection and operation to a specific surgical smoke evacuator, please refer to the specific instructions for use that accompany that system for appropriate set-up and use.

Prior to connecting the PlumePen® Pro 3-prong power cord to the electro-surgical generator, ensure that the electrode blade (provided) is fully inserted into the pencil by gently pressing it in front of the PlumePen® Pro. Remove the plastic shield from the blade and discard.

Figure 2. Plug the 3-prong power cord into the electrosurgical generator of your choice, (Figure 2). Confirm that all power settings on the generator are appropriate for the procedure being performed.

The electrosurgical generator's intensity should be set as low as is necessary to achieve the desired effect. PlumePen® Pro is a monopolar electrode, the use of a dispersive electrode is required to prevent burns/injury to patient.

Please refer to electrosurgical generator user manual and dispersive electrode instructions for use for additional instructions.


Figure 3. PlumePen® Pro is equipped with several key features for the surgeon's convenience. Surrounding the electrode blade is a translucent capture port tip that suctions surgical smoke as it is created.

Keep the active electrodes clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning; injury to operating room personnel may result.



Figure 4. PlumePen® Pro has CUT and COAG buttons located on the top of the electrosurgical pencil body. The YELLOW button is for cutting and the BLUE button is for coagulation, (Figure 4).

After use, turn off surgical smoke evacuator and electrosurgical generator. Discard entire PlumePen® Pro assembly (pencil, blade, tubing, and power cord) per facility protocol for contaminated items.

Figures 5 / 6. Turn off generator. To remove capture port tip, press firmly in a downward motion to release as shown in Figures 5 and 6.


Capture port tip contains a radiopaque marker. Remove capture port tip in a location not directly over an open incision.

Figures 7 / 8. To remove blade, grasp blade at hub and pull forward as shown in Figures 7 and 8.

Figures 9 / 10. **CAUTION:** We recommend the use of the original Buffalo Filter hexed blade or a blade containing a hexagonal hub to prevent blade movement. If original blade is removed, confirm new blade is completely inserted and secured before activating the pencil. Never force the blade into the pencil.

Figures 11 / 12. To replace capture port tip, align bottom tab as shown in Figure 11. Then press down on capture port tip tab to lock in place as shown in Figure 12.

CAUTION:

DO NOT USE in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.

Sterile unless packaging is damaged or any seal is broken. Do not use electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids, or objects, or in the presence of oxidizing agents, as fire may result. When not in use, the active electrode should be placed in a clean, dry, non-conductive safety holster. Inadvertent contact with the patient may result in burns. Contact with drapes or linens may cause a fire. DO NOT activate the instrument when not in contact with target tissue, as this may cause injuries due to capacitive coupling.

PlumePen® Pro is designed and intended only to be used with an electrosurgical generator that has been tested to the IEC 60601 standard. Please refer to generator to ensure compatibility. PlumePen® Pro is not a fluid removal device, therefore, it should not be used for such application.

ASPIRATE fluid from the area before activating the instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient.

For procedures where visualization may be impaired, be alert of these potential hazards:


- The electrode tip may remain hot enough to cause burns after the current has been de-activated.
- Inadvertent activation or movement of the activated electrode outside of the field of vision may result in injury to the patient.
- Localized burns to the patient or physician may result from electrical currents carried through conductive objects. Current may be generated in conductive objects by direct contact with the active electrode, or by the active accessory being in close proximity to the conductive object.

For single use only. Do not reuse, reprocess, or re-sterilize. Reuse, reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness, or death.

If you have any questions or need any assistance with the set-up or use of PlumePen® Pro, please contact Buffalo Filter at:

5900 Genesee Street, Lancaster, NY 14086
 U.S. toll-free: 800.343.2324, Fax: 716.835.3414
 International: 1.716.835.7000; Email: smoke-info@buffalofilter.com
 www.buffalofilter.com
 904316 RevB


Not Made with Natural Rubber Latex

Do Not Use if Package is Damaged

Single Use Only Do Not Reuse

Emergo Europe
 Molenstraat 15
 2513 BH, The Hague
 The Netherlands

REF PLPRO4000
 PLPRO4500

CAUTION:
 Intended for use with a maximum voltage of 10kV P-P.

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