

# VisiClear®

## Surgical Smoke Evacuator



### Procedure Operating Modes

*One of three operating modes (Laparoscopic, Open Tubing, and Electrosurgical Pencil Mode) can easily be selected, automatically customizing filter life based on surgical procedure.*

### Advanced Sound Control Design

*Whisper Technology™ ensures ultra-quiet operation with unprecedented sound quality.*

### Sleek Shelf Design

*Ergonomically designed to accommodate stacking.*

### Occlusion Warning

*Safety feature that alerts medical personnel when an occlusion has occurred.*

### Intuitive Color Touch Screen

*Easy to use, requires minimal training, and can be operated with gloved hand.*

### Unique Sliding Port Closure


*SafePort Technology™ ensures reduced risk of staff coming in contact with a contaminated port.*



# VisiClear®

## Surgical Smoke Evacuator

VisiClear® is our most advanced acute care offering for surgical smoke evacuation. Combining Safeport Technology™, an occlusion warning, filter life tracking, and one-touch programmable procedure modes for open tubing, electrosurgical pencil, and laparoscopic procedures, VisiClear effectively captures and filters surgical smoke, removes odors, particulates and other potentially hazardous by-products.



### Technical Information

**Noise Levels** (Nominal Values)  
 @ minimum power < 43 dBA  
 @ maximum power < 55 dBA

**Dimensions** (H x W x D)  
 7.1 in x 13.9 in x 16.1 in  
 18.0 cm x 35.4 cm x 40.9 cm

**Flow Rates**  
 30 cfm (839 lpm) using 7/8 in (22 mm) tubing

**Weight**  
 With filter installed: 17.5 lbs (7.9 kg)  
 Without filter installed: 14.5 lbs (6.6 kg)

**Static Suction**  
 150 mmHg vacuum

**Certifications**  
 UL, cUI, CE

**Warranty**  
 Three (3) Year Full Warranty

### Product Ordering Information

Description	Catalog Number	Quantity
VisiClear® (100/120 Volt Systems) Includes one filter and one footswitch	VC120	1/Each
VisiClear® (200/240 Volt Systems) Includes one filter and one footswitch	VC220	1/Each
VisiClear® (100/120 Volt Systems) Includes one footswitch	VC120 NF	1/Each
VisiClear® (200/240 Volt Systems) Includes one footswitch	VC220 NF	1/Each
EZLink® Automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life.	EZLINK01	1/Case
ViroSafe® Fluid Trap Designed to capture fluids to help maintain the integrity of the filtration system and assist with maintaining filter life.	VSFT10	10/Case
ViroSafe® VisiClear® Filter	VS135	1/Case
Safeport Technology™ prevents the clinician from accidentally touching an exposed port when not in use.	VS1352	2/Case
	VS1354	4/Case

