

Surgical Smoke Plume Evacuators	VisiClear®	ViroVac®	PlumeSafe® Turbo	Porta PlumeSafe® 604
Market Segment	Acute Care	Outpatient Surgery Centers & Physicians Offices	Aesthetic	Non-acute Care & Veterinary
Procedure Operating Modes: - Laparoscopic - Electrosurgical Pencil - Open Tubing	Yes	No	No	No
Sound / Noise Levels (decibels @ 25 CFM)	51 dBA	55 dBA	50 dBA	Not available
Stackability / Shelf Friendly	Yes	Yes	No	No
Occlusion Warning	Yes	No	No	No
Display Type	Color Touch Screen	LED Display	LCD Display	LED Display
3-port Filter - 1/4 in (6.4 mm) - 3/8 in (9.5 mm) - 7/8 in (22 mm)	Yes	Yes	No	No
CE Marked	Yes	Yes	Yes	No
ROHS Compliant	Yes	Yes	Yes	No