



Buffalo Filter, LLC
 5900 Genesee Street, Lancaster, NY 14086, U.S.A.
 Phone: 800.343.2324 (US); +1.716.835.7000; Fax: 716.835.3414
 Email: bf.smoke.info@buffalofilter.com Web: www.buffalofilter.com

PlumePort® Laparoscopic Devices	PlumePort® ActiV®	PlumePort® SEO	PlumePort® SE	PlumePort® SE	PlumePort®
Type of Device	Active	Active	Active	Active	Passive
Suction Method	Standard suction	Smoke evacuator	Smoke evacuator	Smoke evacuator	Not applicable
Flow Adjustment	Yes	Via smoke plume evacuator	Via smoke plume evacuator	Via smoke plume evacuator	Yes
Flow Adjustment Mechanism	Flow Adjustment Dial on sterile field	Calibrated orifice requires no OR personnel interaction	Controlled by trocar stopcock (not recommended with BF ViroVac)	Controlled by valve on sterile tubing (limits CO2 usage)	Flow Adjustment Dial on sterile field
Sound	Noiseless	Based on smoke evacuator flow	Based on smoke evacuator flow	Based on smoke evacuator flow	Noiseless
Built in Fluid Trap	Yes	No	No	No	Yes
Latex Free	Yes	Yes	Yes	Yes	Yes
Non-Phthalate Tubing	Yes	No	No	No	No