

# Patents

Product Description	Catalog Number	Patent Number	Quantity	Sterile/ Non-Sterile
<b>Surgical Smoke Evacuators</b>				
<b>VisiClear® Surgical Smoke Evacuator</b>				
100/120 Volt Systems	VC120	7,761,188 8,095,241	-	-
220/240 Volt Systems (includes one filter and one footswitch)	VC220		-	-
100/120 Volt Systems	VC120NF		-	-
220/240 Volt Systems (without filter, includes one footswitch)	VC220NF		-	-
ViroSafe® VisiClear® Filter	VS135 Series		-	-
ViroVac® Surgical Smoke Evacuator			8,147,577 7,761,188 8,095,241	-
100/120 Volt Systems	VV120	-		-
220/240 Volt Systems (includes one filter and one footswitch)	VV220	-		-
100/120 Volt Systems	VV120NF	-		-
220/240 Volt Systems (without filter, includes one footswitch)	VV220NF	-		-
ViroSafe® ViroVac® Filter	VS353 Series	-		-
<b>PlumeSafe® Turbo Surgical Smoke Evacuator</b>				
100/120 Volt Systems	PSWTURBO AS	7,761,188 8,095,241	-	-
220/240 Volt Systems (includes one footswitch; filter sold separately)	PSWTURBO AS/220		-	-
PlumeSafe® Turbo Replacement Filters	VS06, VSXLL, VSXLLT, VSEY, VSOH, & VSVET Series		-	-
<b>Porta PlumeSafe® 604 Surgical Smoke Evacuator</b>				
100/120 Volt Systems (includes one footswitch; filter sold separately)	PPS604	7,761,188 8,095,241	-	-
Porta PlumeSafe® 604 Replacement Filters	VS06, VSXLL, VSXLLT, VSEY, VSOH, & VSVET Series		-	-
<b>ViroSafe® Filters</b>				
<b>ViroSafe® 6 Filter</b>				
	VS06001	7,761,188 8,095,241	1/Case	-
	VS06002		2/Case	-
	VS06004		4/Case	-
	VS06012		12/Case	-
<b>ViroSafe® 6 Twin Port Filter</b>	VS06T04	7,761,188 8,095,241	4/Case	-
<b>ViroSafe® Extra Long Life Filter</b>				
	VSXLL02	7,761,188 8,095,241	2/Case	-
	VSXLL04		4/Case	-
<b>ViroSafe® Extra Long Life Twin Port Filter</b>				
	SINGVSXLLT	7,761,188 8,095,241	1/Case	-
	VSXLLT02		2/Case	-
	VSXLLT04		4/Case	-
<b>ViroSafe® Veterinary Filter</b>	VSVET01	7,761,188 8,095,241	1/Case	-
<b>ViroSafe® Erbium: YAG Pre-filter</b>	BFEY04PF	-	4/Case	-
<b>ViroSafe® Erbium: YAG Filter</b>	SING VSEY	7,761,188 8,095,241	1/Case	-

# Patents

Product Description	Catalog Number	Patent Number	Quantity	Sterile/ Non-Sterile
ViroSafe® Erbium: YAG Kit	VSEY1-4	7,761,188 8,095,241	1/Main Filter 4/Pre-filters	-
ViroSafe® Odor/Hair Removal Filter	VSOH001	7,761,188 8,095,241	1/Case	-
	VSOH002		2/Case	-
ViroSafe® Pre-Filter	BFPF7/8-10	-	10/Case	-
ViroSafe® 12 Filter	VS12D01	-	1/Case	-
	VS12D02	-	2/Case	-
ViroSafe® ViroVac® Filter	VS353	8,147,577 7,761,188 8,095,241	1/Case	-
	VS35302		2/Case	-
	VS35304		4/Case	-
ViroSafe® VisiClear® Filter	VS135	7,761,188 8,095,241	1/Case	-
	VS1352		2/Case	-
	VS1354		4/Case	-
<b>In-Line Filters</b>				
BILF150 Multi-Use In-Line Filter	BILF150	-	10/Case	-
BILF200 Multi-Use In-Line Filter	BILF200	-	10/Case	-
BILF100 Single Use In-Line Filter	BILF100	-	10/Case	-
BILF101 Single Use In-Line Filter	BILF101	-	10/Case	-
<b>Competitive Replacement Filters</b>				<b>Original Filter</b>
Berchtold Next Generation Phoenix XS, X3	BTLDXS04	7,597,731 7,892,337	4/Case	Berchtold Phoenix XS, X3
Berchtold TELEVAC Systems	BTLD004	7,597,731 7,892,337	4/Case	Berchtold TELEVAC, TELEVAC LT; Surgimedics 901428.
Cooper Surgical System 6080	BSUCF03V	-	3/Case	Coastal Life 1101300, 901300-000; Cooper 6082; Intermed IG-6082; Surgimedics 901300, 901009; Valleylab E3515.
Cooper Surgical System 6080; Surgimedics: SurgiFresh, Tabletop, Mini, & Plume-inator Systems; Valleylab AirForce System	BSUCF03	-	3/Case	Coastal Life 1101300, 901300-000; Cooper 6082; Intermed IG-6082; Surgimedics 901300, 901009; Valleylab E3515.
Cooper Surgical System 6080; Surgimedics: SurgiFresh, Tabletop, Mini, & Plume-inator Systems; Valleylab AirForce System	BSPF7-8-30	-	30/Case	Cooper 6081; Surgimedics 901202; Valleylab E3510, E3630.
Cooper Surgical System 6080; Surgimedics: SurgiFresh, Tabletop, Mini, & Plume-inator Systems; Valleylab: AirForce, Optimumm Systems	BSPF030	-	30/Case	Coastal Life 1101202, 901202-000; Cooper 6081; Intermed 2081, 6081; Surgimedics 901202; Valleylab E3510, E3630.
EDGE S.A.F.E. System	BFE20004	-	5/Case	EDGE 20004

# Patents

Product Description	Catalog Number	Patent Number	Quantity	Sterile/ Non-Sterile
EDGE S.A.F.E. System	BFE20009	-	1/Case	EDGE 20009
EDGE S.A.F.E. System	BFE Kit	-	1/1 main filter; 5 pre-filters /Case	EDGE 20009, 20004
QuantaVac System	MFSD8-1	7,761,188 8,095,241	1/Case	QuantaVac System main filter.
Stackhouse AirSafe VersaVac I & II	BST370ULC	-	3/Case	EDGE 20012-03; Stackhouse ST-370-FRESH/3 (24300); Surgimedics 901408.
Stackhouse ES-2000 System	BES2000ULC	-	1/Case	Stackhouse ES-2000/ULC (24152); Surgimedics 901407.
Stackhouse ES-2000 System	BES2000PFC	-	12/Case	Lumenis 8000-002; Stackhouse ES-2000PF/12 (24204), ES-20000PFC/12 (24205); Surgimedics 901405, 901406.
Stackhouse ES-2000 System	BESIWT110	-	12/Case	Stackhouse ES-IWT-110/S (24211); Surgimedics 901414.
Stackhouse ES-2000 System	BES2000PF	-	12/Case	Lumenis 8000-002; Stackhouse ES-2000PF/12 (24204), ES-20000PFC/12 (24205); Surgimedics 901405, 901406.
Stackhouse MiniVac ST-1250	BSTM300	-	3/Case	Conmed 60-5902-001, Cooper 395-703, Stackhouse STM-300 (24139-01); Surgimedics 901409; Valleylab E1091.
Stackhouse: LFA-100, Point One Systems	BPTU100	-	1/Case	Stackhouse PTU-100 (24000); Surgimedics 901411.
Stackhouse: LFA-100, Point One Systems	BLFS103	-	3/Case	Stackhouse LFS-103/3 (24001); Surgimedics 901410.
Stackhouse: LFA-100, Point One Systems	BLFS101	-	12/Case	Stackhouse LFS-103/12 (24033); Surgimedics 901410.
Stryker Neptune 2 System	VSNEP01	8,147,577	1/Case	Stryker Neptune 2 Ultra
Stryker Neptune 2 System	VSNEP04	8,147,577	4/Case	Stryker Neptune 2 Ultra
StrykeVac 1 & 2 Smoke Evacuation System	VSSV505	8,147,577 7,761,188 8,095,241 7,597,731 7,892,337	1/Case	StrykeVac ViroSafe
StrykeVac 1 & 2 Smoke Evacuation System	VSSV50502	8,147,577 7,761,188 8,095,241 7,597,731 7,892,337	2/Case	StrykeVac ViroSafe
StrykeVac 3 Smoke Evacuation System	VSSV301	8,147,577 7,761,188 8,095,241	1/Case	StrykeVac 3 ViroSafe

# Patents

Product Description	Catalog Number	Patent Number	Quantity	Sterile/ Non-Sterile
StrykeVac 3 Smoke Evacuation System	VSSV302	8,147,577 7,761,188 8,095,241	2/Case	StrykeVac 3 ViroSafe
StrykeVac 3 Smoke Evacuation System	VSSV304	8,147,577 7,761,188 8,095,241	4/Case	StrykeVac 3 ViroSafe
TRUMPF EasyVac System	VSTR08	8,147,577 7,761,188 8,095,241 7,597,731 7,892,337	1/Case	TRUMPF ViroSafe Filter
TRUMPF EasyVac System	VSTR0804	8,147,577 7,761,188 8,095,241 7,597,731 7,892,337	4/Case	TRUMPF ViroSafe Filter
Valleylab Optimumm System	BFE3630	-	25/Case	Surgimedics 901202; Valleylab E3630
Valleylab Optimumm System	BFE3625	-	3/Case	Valleylab Optimumm E3625; Surgimedics 901412.
<b>Surgical Smoke Evacuation Pencils and Adapters</b>				
<b>PlumePen® Elite (Stainless Steel ESU Blade)</b>	PLP1005	9,289,261	5/Case	-
Includes 3/8 in (9.5 mm) x 10 ft (3 m) tubing with adapter to 7/8 in (22 mm) port.	PLP1020		20/Case	-
<b>PlumePen® Elite (Stainless Steel ESU Blade)</b>	PLP1505	9,289,261	5/Case	-
Includes 3/8 in (9.5 mm) x 15 ft (4.6 m) tubing with adapter to 7/8 in (22 mm) port.	PLP1520		20/Case	-
<b>PlumePen® Elite (Non-Stick ESU Blade)</b>	PLP2005	9,289,261	5/Case	-
Includes 3/8 in (9.5 mm) x 10 ft (3 m) tubing with adapter to 7/8 in (22 mm) port.	PLP2020		20/Case	-
<b>PlumePen® Elite (Non-Stick ESU Blade)</b>	PLP2505	9,289,261	5/Case	-
Includes 3/8 in (9.5 mm) x 15 ft (4.6 m) tubing with adapter to 7/8 in (22 mm) port.	PLP2520		20/Case	-
<b>PlumePen® Pro (Non-Stick ESU Blade)</b>	PLPRO4020	9,289,261	20/Case	-
Includes 3/8 in (9.5 mm) x 10 ft (3 m) tubing with adapter to 7/8 in (22 mm) port.				-
<b>PlumePen® Pro (Non-Stick ESU Blade)</b>	PLPRO4520	9,289,261	20/Case	-
Includes 3/8 in (9.5 mm) x 15 ft (4.6 m) tubing with adapter to 7/8 in (22 mm) port.				-
<b>PlumePen® Pro Capture Ports (x-ray detectable):</b> (Accommodates 4 in [102 mm] electrodes)	PROTIP420	9,289,261	20/Case	-
<b>PlumePen® Pro Capture Ports (x-ray detectable):</b> (Accommodates 6 in [152 mm] electrodes)	PROTIP620	9,289,261	20/Case	-
<b>PenAdapt®</b>	PA1025	6,524,307	25/Case	-
Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.				
<b>PenAdapt® (Cut-back)</b>	PA1025CB	6,524,307	25/Case	-
Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.				
<b>PenAdapt®</b>	PA2010	6,524,307	10/Case	-
Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.83 m) tube.				
<b>PenAdapt® (Cut-back)</b>	PA2010CB	6,524,307	10/Case	-
Cut-back adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.83 m) tube.				
<b>PenAdapt®</b>	PA2010B	6,524,307	10/Case	-
Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.83 m) tube.				
<b>Laparoscopic Solutions</b>				
<b>PlumePort® ActiV® Laparoscopic Smoke Filtration Device</b>	PP010CS	8,608,816 9,415,160	10/Case	-
<b>PlumePort® SEO Laparoscopic Smoke Filtration Device</b>	PP012SEO	-	12/Case	-

# Patents

Product Description	Catalog Number	Patent Number	Quantity	Sterile/ Non-Sterile
<b>PlumePort® SEV Laparoscopic Smoke Filtration Device</b> Includes 7/8 in (22 mm) filter inlet, 1/4 in (6 mm) x 10 ft (3 m) with surgical smoke control valve.	BLAP100V	-	12/Case	-
<b>PlumePort® SE Laparoscopic Smoke Filtration Device</b> Includes 7/8 in (22 mm) filter inlet, 1/4 in (6 mm) x 10 ft (3 m).	BLAP100	-	12/Case	-
<b>PlumePort® Laparoscopic Smoke Filtration Device</b>	PP010	8,608,816 9,289,261 9,011,366	10/Case	-
<b>LaparoVue® Visibility System</b>	LAPVUE10	-	10/Case	-
<b>LaparoLight™ Veress Needle</b>	LLVN12012 (120mm) LLVN15012 (150mm)	8,523,817 8,838,206	12/Case	-
<b>Tubing - Standard &amp; Specialty</b>				
1/4 in (6.4 mm) x 2 ft (.61 m).	ST10105	-	5/Case	Non-Sterile
3/8 in (9.5 mm) x 10 ft (3.1 m).	VT38S1010	-	10/Case	Sterile
7/8 in (22 mm) x 6 ft (1.8 m).	VT10424	-	24/Case	Sterile
7/8 in (22 mm) x 10 ft (3.1 m).	VT10624	-	24/Case	Sterile
7/8 in (22 mm) x 6 ft (1.8 m) with (12) 7/8 in tubing-to-tubing connectors (extension tubing).	VT10324	-	24/Case	Non-Sterile
7/8 in (22 mm) x 10 ft (3.1 m).	VT10524-10	-	10/Case	Non-Sterile
7/8 in (22 mm) x 10 ft (3.1 m).	VT10524	-	24/Case	Non-Sterile
7/8 in (22 mm) x 6 ft (1.8 m).	VT78-6	-	6/Case	Non-Sterile
7/8 in (22 mm) x 2.5 ft (.76 m) (Replaces suction canister connecting tube).	VT10610	-	10/Case	Non-Sterile
1-1/4 in (32 mm) x 6 ft (1.8 m).	VT10206	-	6/Case	Sterile
1-1/4 in (32 mm) x 6 ft (1.8 m).	VT10106	-	6/Case	Non-Sterile
1-1/4 in (32 mm) x 10 ft (3.1 m).	VT10306	-	6/Case	Non-Sterile
1-1/4 in (32mm) x 2 ft (.61 m).	VT10510	-	10/Case	Non-Sterile
1-3/8 in (35 mm) x 6 ft (1.85 m) tube (replace as needed) and six 2 ft (.61 m) length tubes with integral wand & sponge guard to be replaced each procedure.	VTEYKIT	-	1/Kit	Non-Sterile
7/8 in (22 mm) x 6 ft (1.8 m) with 7/8 in (22 mm) to 1/4 in (6.4 mm) reducer to 1/4 in (6.4 mm) x 2 ft (.6 m) tubing.	VTESU05	-	5/Case	Sterile
7/8 in (22mm) x 6 ft (1.8 m) with 7/8 in (22 mm) to 1/4 in (6.4 mm) reducer to 1/4 in (6.4 cm) x 2 ft (.6 m) tubing.	PSFTS05	-	5/Case	Non-Sterile
7/8 in (22 mm) x 6 ft (1.8 m) with 7/8 in (22 mm) to 3/8 in (9.5 mm) reducer and 3/8 in (9.5 mm) x 4 ft (1.2 m).	VTVIC10	-	10/Case	Sterile
7/8 in (22 mm) to 3/8 in (9.5 mm) reducer fitting with 7/8 in (22 cm) x 6 ft (1.8 m) tubing.	VT10405	-	5/Case	Sterile
7/8 in (22 mm) to 1/4 in (6.4 mm) reducer fitting with x 2 ft (.6 m) 1/4 in (6.4 mm) tubing.	RF10210	-	10/Case	Sterile
7/8 in (22 mm) to 1/4 in (6.4 mm) reducer fitting with x 2 ft (.6 m) 1/4 in (6.4 mm) tubing.	RF10110	-	10/Case	Non-Sterile
7/8 in (22 mm) x 6 ft (1.8 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT424	-	24/Case	Sterile
7/8 in (22 mm) x 10 ft (3.1 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT624	-	24/Case	Sterile
7/8 in (22 mm) x 14 ft (4.3 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT805	-	5/Case	Sterile
7/8 in (22 mm) x 10 ft (3.1 m) tubing with 8 in (20 cm) integral wand and sponge guard with 7/8 in to 1-3/8 in (3.5 cm) adapter.	VTWTS10	-	10/Case	Sterile

# Patents

Product Description	Catalog Number	Patent Number	Quantity	Sterile/ Non-Sterile
7/8 in (22 mm) x 6 ft (1.8 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT324	-	24/Case	Non-Sterile
7/8 in (22 mm) x 10 ft (3.1 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT524	-	24/Case	Non-Sterile
1-1/4 in (32 mm) x 6 ft (1.8 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT206	-	6/Case	Sterile
1-1/4 in (32 mm) x 10 ft (3.1 m) tubing with 8 in (20 cm) integral wand and sponge guard.	VTWT406	-	6/Case	Sterile
1-3/8 in (35 mm) x 6 ft (1.8 m) with wand and sponge guard.	VTWT505	-	5/Case	Sterile
1-3/8 in (35 mm) x 10 ft (3.1 m) with wand and sponge guard.	VTWT515	-	5/Case	Sterile
1-3/8 in (35 mm) x 6 ft (1.8 m) with wand and sponge guard.	VTWT506	-	6/Case	Non-Sterile
<b>Tubing - Reducer Fittings/Adapters</b>				
7/8 in (22 mm) to 1/4 in (6.4 mm) reducer fitting.	RF10112	-	12/Case	Non-Sterile
7/8 in (22 mm) to 1/4 in (6.4 mm) reducer fitting with 2 in (5.1 cm) cuff.	RF12524	-	24/Case	Non-Sterile
7/8 in. (22 mm) to 3/8 in (9.5 mm) reducer fitting with 2 in (5.1 cm) cuff.	RF783824	-	24/Case	Non-Sterile
7/8 in (22 mm) to 3/8 in (9.5 cm) reducer fitting.	RF10512	-	12/Case	Non-Sterile
1-3/8 in (35 mm) to 7/8 in (22 mm) reducer fitting.	RF10705	-	5/Case	Non-Sterile
1-3/8 in (35 mm) to 1/4 in (6.4 mm) reducer fitting.	RF10905	-	5/Case	Non-Sterile
1-3/8 in (35 mm) to 3/8 in (9.5 mm) reducer fitting.	RF13810	-	10/Case	Non-Sterile
<b>Tubing - Wands &amp; Sponge Guards</b>				
1-1/4 in (32 mm) x 8 in (20 cm).	WLRT103	-	3/Case	Sterile
7/8 in (22 mm) x 8 in (20 cm).	WLRT410	-	10/Case	Sterile
1-1/4 in (32 mm) x 8 in (20 cm).	WLRT003	-	3/Case	Non-Sterile
7/8 in (22 mm) x 8 in (20 cm).	WLRT303	-	3/Case	Non-Sterile
Sponge guard fits 7/8 in (22 mm) tubing.	WTIPS10	-	10/Case	-
Sponge guard fits 1-1/4 in (32 mm) to 1-3/8 in (35 mm) tubing.	WTIPL10	-	10/Case	-
<b>Accessories</b>				
Armstand	ARMSTAND	-	1/Each	-
Canister, Fluid Suction Collection	SCC1001	-	6/Case	-
EvacuCart™	ETC01	-	1/Each	-
EZLink®	EZLINK01	9,164,500	1/Case	-
LaseAdapt® ProFractional Therapy Attachment	LA3010PF10	8,641,488	10/Case	-
LaseAdapt® SC Contour TRL Erbium Scanner Attachment	LA3010EY10	8,641,488	10/Case	-
LaseAdapt® Reuseable Tubing Clips	LA3010CL10	-	4/Case	-
PlumePen® Pro Capture Ports (x-ray detectable)	PROTIP420	-	20/Case	-
	PROTIP620	-	20/Case	-
Pneumatic Footswitch	FS10004	-	1/Case	-
Ring Holder for Suction Canister	SCR1001	-	1/Case	-
Surgical Laser Masks	SMLOF50	-	50/Case	-
	SMLOF300	-	300/Case	-
ViroSafe® Fluid Trap	VSFT10	-	10/Case	-
ViroVac® Wall Mount Kit	VVWM01	-	1/Each	-
VueTip™ Trocar Swab Accessory Kit	VueTip01	-	1/Each	-
	VueTip10	-	10/Case	-
Wall Plate for Suction Canister Ring	WALLPLATE	-	1/Case	-
Y Adapter	ADP22X22Y	-	1/Each	-