



LaparoVue® Visibility System



*change your Vue...
with an All-in-One Solution*



Product Ordering Information

Description	Catalog Number	Quantity
LaparoVue® Visibility System (includes main unit/scope cradle, 2 trocar swabs, and 1 microfiber cloth)	LAPVUE10	10/Case
VueTip™ Trocar Swab Accessory Kit (includes 2 trocar swabs and 1 microfiber cloth)	VUETIP10	10/Case



VueTip™ Trocar Swab Accessory Kit
(sold separately; VUETIP10)



Trial Now — call 1.716.835.7000 or 1.800.343.2324

To view additional **Laparoscopic Solutions** visit www.buffalofilter.com.

LaparoVue® Visibility System

Laparoscopic Solutions



All-in-One Solution
for improved laparoscopic visibility

Laparoscopic Visualization

Complexity Made Simple

Distinctive By Design

warms & white balances

- warms scope in 5 minutes to prevent fogging
- unit warms up to 5 hours
- allows for white balancing



Unique Product Attributes

one unit, two ports

- two ports: warming / white balancing and cleaning / defogging
- fits 3mm–12mm scopes without port size converter
- allows for simultaneous use of both ports

detachable scope cradle

- safely and securely supports scope
- reduces tension on lens and potential for damage

built-in LED lights

- provides enhanced port visibility

ergonomic design

- stable base prevents tipping
- contoured, easy grip profile
- engineered with soft silicone material



cleans

- provides crisp, clear view in 5 seconds
- minimizes procedural delays
- radiopaque VueTip™ Trocar Swabs remove debris and excess liquid from trocar cannula
- x-ray detectable Microfiber Cloth designed to capture residue and particulates at a microscopic level

defogs

- pre-filled warming solution helps prevent intra-operative scope fogging



* Photographs are representative of a fogged and clear / defogged lens but do not reflect actual product use.

Actual product may differ slightly in appearance to images shown.