



Buffalo Filter, LLC  
 5900 Genesee Street, Lancaster, NY 14086, U.S.A.  
 Phone: 800.343.2324 (US); +1.716.835.7000; Fax: 716.835.3414  
 Email: bf.smoke.info@buffalofilter.com Web: www.buffalofilter.com

Surgical Smoke Plume Evacuators	VisiClear®	ViroVac®	PlumeSafe® Turbo	Porta PlumeSafe® 604
<b>Market Segment</b>	Acute Care	Outpatient Surgery Centers & Physicians Offices	Aesthetic	Non-acute Care & Veterinary
<b>Procedure Operating Modes:</b> - Laparoscopic - Electrosurgical Pencil - Open Tubing	Yes	No	No	No
<b>Sound / Noise Levels (decibels @ 25 CFM)</b>	51 dBA	55 dBA	50 dBA	Not available
<b>Stackability / Shelf Friendly</b>	Yes	Yes	No	No
<b>Occlusion Warning</b>	Yes	No	No	No
<b>Display Type</b>	Color Touch Screen	LED Display	LCD Display	LED Display
<b>3-port Filter</b> - 1/4 in (6.4 mm) - 3/8 in (9.5 mm) - 7/8 in (22 mm)	Yes	Yes	No	No
<b>CE Marked</b>	Yes	Yes	Yes	No
<b>ROHS Compliant</b>	Yes	Yes	Yes	No